

Simple Install Guide

Page 1 / 4

Smart Prosthetics ERP · Local network · No internet needed after setup

In simple words: One computer = **server** (holds the program + data). Other computers only open **Chrome** and use a link. No program install on staff PCs.

Client name: _____

Install date: _____

Your name: _____

Server IP (write here):

Staff link: http://_____
/prosthetics/public

DB password: _____

Part A — At home (with internet) — prepare USB

1. Download **Laragon Full** → laragon.org/download
2. Download **Git** (optional) → git-scm.com/download/win
3. On your PC with Laragon, run:

```
cd C:\laragon\www
git clone https://github.com/Refaat1942/pros.git prosthetics
cd prosthetics
git pull origin master
composer install --no-dev --optimize-autoloader
```

PHP ext-zip required: In Laragon → PHP → php.ini, enable extension=zip, then restart Laragon before running Composer. Excel catalog import/export needs this extension.

Important: Copy the whole folder C:\laragon\www\prosthetics to USB. It must include the **vendor** folder. Without it, the system will not work offline.

Part B — At the client — the server PC

Step 1 — Choose the server computer

- One Windows PC that stays ON during work hours
- Connect a UPS (battery backup) if possible
- Ask IT for a fixed IP (example: 192.168.1.50)

Step 2 — Install Laragon

1. Run Laragon installer from USB → install to C:\laragon
2. Open Laragon → click **Start All** (Apache must be green)

Step 3 — Install PostgreSQL (database)

1. Laragon → **Menu** → **PostgreSQL** → **Install**
2. Click **Start All** again

PostgreSQL = where all patient and system data is stored.

Copy project & create database

Page 2 / 4

Steps 4–7

Step 4 — Copy the project from USB

Copy folder to: C:\laragon\www\prosthetics

Check inside: folders vendor, public, file .env.example

Step 5 — Create database (one time)

Laragon → **Menu** → **PostgreSQL** → **psql** (black window opens)

Paste this (change the password!):

```
CREATE DATABASE prosthetics ENCODING 'UTF8' LC_COLLATE 'C' LC_CTYPE 'C' TEMPLATE
template0;
CREATE USER prosthetics_user WITH PASSWORD 'StrongPassword2026!ChangeMe';
GRANT ALL PRIVILEGES ON DATABASE prosthetics TO prosthetics_user;
\c prosthetics
GRANT ALL ON SCHEMA public TO prosthetics_user;
GRANT CREATE ON SCHEMA public TO prosthetics_user;
\q
```

\q = exit. Write the same password in Step 6.

Step 6 — Edit settings file (.env)

Open Command Prompt (search "cmd") and type:

```
cd C:\laragon\www\prosthetics
copy .env.example .env
php artisan key:generate
notepad .env
```

In Notepad, change these lines (use real server IP):

```
APP_URL=http://192.168.1.50
APP_DEBUG=false
DB_DATABASE=prosthetics
DB_USERNAME=prosthetics_user
DB_PASSWORD=StrongPassword2026!ChangeMe
SESSION_SECURE_COOKIE=false
TELEGRAM_ERROR_NOTIFY=false
FIREBASE_ENABLED=false
OFFLINE_ASSETS=true
```

Save and close Notepad.

Step 7 — Start the system (run once)

```
cd C:\laragon\www\prosthetics
php artisan migrate --force
php artisan db:seed --class=Database\Seeders\RolesAndAdminSeeder
php artisan storage:link
php artisan config:cache
```

Command	What it does (simple)
migrate	Creates empty tables in the database
db:seed RolesAndAdminSeeder	Creates login accounts (admin, reception, doctor...)

storage:link	Lets system save uploaded files
config:cache	Saves settings for faster running

Do **NOT** run full `php artisan db:seed` — it adds fake demo patients. Use **RolesAndAdminSeeder** only.

Login, network & staff PCs

Page 3 / 4

Steps 8–11

Step 8 — First login (on server)

1. Open **Chrome**
2. Go to: `http://prosthetics.test` OR `http://YOUR-IP/prosthetics/public`
3. Login: username **superadmin** · password **123456**
4. **Change password immediately!**

Step 9 — Test from another computer

On another PC on the same network, open Chrome:

```
http://192.168.1.50/prosthetics/public
```

(Replace with your real server IP)

If it does not work:

- Server is ON and Laragon → Start All (green)
- Windows Firewall: allow Apache on Private network (ask IT)
- Wrong IP — run `ipconfig` on server, find IPv4 Address

Step 10 — Staff computers (easy)

On each staff PC — **nothing to install**:

1. Open Chrome → go to the link above
2. Bookmark (Ctrl+D)
3. Optional: desktop shortcut named "Prosthetics ERP"

Step 11 — Backup (important)

Test once on server:

```
cd C:\laragon\www\prosthetics
php artisan prosthetics:backup
```

Backup files: `C:\laragon\www\prosthetics\storage\backups`

Schedule daily at 2:00 AM (Task Scheduler) or copy folder to USB weekly.

Words you can ignore (simple meaning)

Word	Meaning
Laravel	The program framework — Laragon runs it for you
PostgreSQL	Database — stores all records
.env	Settings file — edit once
migrate	Build empty tables in database
vendor	Pre-downloaded libraries — must be on USB
Composer	Downloads libraries — do at home, not at client

Checklists & troubleshooting

Page 4 / 4

Print and use on site

Before visit checklist

- ☐ USB has Laragon installer
- ☐ USB has full prosthetics folder (with vendor)
- ☐ Wrote client name, date, your phone number

On server checklist

- ☐ Laragon installed → Start All green
- ☐ PostgreSQL installed → Start All green
- ☐ Project in C:\laragon\www\prosthetics
- ☐ Database created (Step 5)
- ☐ .env edited (IP + password)
- ☐ Step 7 commands finished OK
- ☐ Login works on server
- ☐ Login works from another PC
- ☐ superadmin password changed

Staff handover checklist

- ☐ Gave everyone the link
- ☐ Created real user accounts (not default 123456)
- ☐ Backup tested
- ☐ Client knows: do not turn off server during work hours

If something breaks

Problem	Try this
Blank page / error	Laragon → Stop All → Start All
Cannot open from another PC	Check IP, Firewall, server is ON
Database error (red text)	Password in .env must match Step 5
After power cut	Turn server ON → Laragon → Start All

Later update (when internet is back — server only)

```
cd C:\laragon\www\prosthetics
git pull origin master
```

```
composer install --no-dev --optimize-autoloader  
php artisan migrate --force  
php artisan config:cache
```

Tell staff: press **Ctrl + F5** in Chrome.

One-line summary: One PC runs Laragon + the program. Everyone else opens a link in Chrome on the local network. No internet needed for daily work.